SMALL VOLUME PACKAGE

SMALL VOLUME PACKAGE SVP 180

(Ref. 600950)

SMALL VOLUME PACKAGE SVP 65

(Ref. 600900)

Temperature range: -20°C to +180°C (with circulation bath not included)



Centring tool ST-R

One spindle MK-SV

One chamber MB-SV

Water jacket

Temperature range: -20°C to +65°C (with circulation bath not included)







- Centring tool ST-R (Ref. 114436).

Tel: +33 (0)4 78 08 54 06 Fax: +33 (0)4 78 08 69 44 contact@lamyrheology.com

11 A, rue des Aulnes 69410 Champagne au Mont d'Or France



www.lamyrheology.com