

# WATER JACKET + SMALL VOLUME PACKAGE

## SMALL VOLUME PACKAGE

### SVP180

(Ref. 600950)



CT-DIN (Ref. 111914), temperature range from -20°C to +180°C.

## SMALL VOLUME PACKAGE

### SVP65

(Ref. 600900)



CT-LC (Ref. 600800), temperature range from -20°C to +65°C.

## WATER JACKET

Our fluid circulation jacket allow use of all measurement systems including cylinder and cup by ensuring precise temperature control associated with a circulation bath (not supplied).

Sold separately or as package, these accessories are compatible with all viscometers and rheometers in the LAMY RHEOLOGY range excepted FIRST PRODIG CP 1000, RM 100 CP 1000/2000 PLUS, RM 200 CP 4000 PLUS, GT-300 PLUS, GT-300 PRODIG and DSR 500 CP 4000 PLUS.

## SMALL VOLUME PACKAGE

### SVP65 OR SVP180

(Ref. 600900 or 600950)

- 1 Water jacket CT-DIN or CT-LC.
- 1 MK-SV spindle choice (p79).
- 1 MB-SV measuring cup choice (p79).
- 1 Centring tool ST-R (Ref. 114436).
- 2 Silicon pipe (1m).

## COMPATIBLES

### MEASURING SYSTEMS

MS-DIN (p72), MS-R according to spindle (p74), MS-VANES according to spindle (p77), MS-ULV (p78), MS-SV (p79).

## DIMENSIONS AND WEIGHT

CT-DIN (Ref. 111914)

Ø100 x H200 mm

Weight 3Kg

CT-LC (Ref. 600800)

D70 x W70 x H160 mm

Weight 1,5Kg

Tel : +33 (0)4 78 08 54 06  
Fax : +33 (0)4 78 08 69 44  
contact@lamyrheology.com

11 A, rue des Aulnes  
69410 Champagne au Mont d'Or  
France



www.lamyrheology.com