

EVA LR PLUS



FEATURES AND BENEFITS

- Temperature range +15°C to +45°C.
- Compatible with bayonet coupling device.
- Temperature setting of 600 ml beakers according to ISO 2555.



600 ml BEAKER



AVAILABLE TEMPERATURE DEVICE

Compatible with bayonet coupling device

Designation	Reference	Temperature range	Type	Compatible measuring system	Dimensions (p x l x h) in mm	Weight
EVA LR PLUS*	N950006	+15°C to +45°C	Peltier	MS-RV/LV	610x340x650	15kg

* Usable with 600ml beaker.